

BE IN THE MOMENT

WITH HANSATON FOKUS



**Rechargeable BTE and
battery-powered RIC
hearing aids**

A Sonova brand

 **HANSATON**
hearing & emotions

HANSATON FOKUS

Great technology, stylish designs




HANSATON FOKUS hearing aids come with a modern, unique design and the exceptional FOKUS sound performance. They're available in five performance levels powered by the innovative HearIntelligence™ technology. They also offer connectivity to all mobile phones* and hands-free phone calls.



HANSATON FOKUS

BATTERY-POWERED RIC HEARING SYSTEM

Battery-powered hearing wonder



At only 2.4 cm, the **sound FS R312** is almost imperceptible and yet impresses with its innovative, powerful technology and a traditional 312 battery.



"These are so tiny and they provide great performance."

CLIENT BENEFIT



Sterling Silver (U4)



Metallic Anthracite (UC)



Carbon Black (U6)







Cashmere Beige (UA)



Mocha Brown (UB)

— Open Dome / Cap Dome
... Vented Dome
--- Power Dome / SlimTip




Receiver type	S	M	P	UP
Output / Gain (2cc)	111 / 47	114 / 51	122 / 59	132 / 71

* Smartphones and traditional mobile phones with compatible Bluetooth® profiles. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.

HANSATON FOKUS

RECHARGEABLE BTE HEARING SYSTEM

Robust, rechargeable, and connected



At 3.3 cm the **AQ beat FS R** impresses with a dual receiver, a rechargeable Lithium-ion battery, and includes tap control for easy access to connectivity.



"The great sound performance and connectivity you want, in a BTE design you will love. Plus easy handling, fast charging, and less waste."

CLIENT BENEFIT



Sterling Silver (U4)



Metallic Anthracite (UC)



Carbon Black (U6)




Cashmere Beige (UA)



Mocha Brown (UB)

Receiver type	P
Output / Gain (2cc) Earhook	130 / 68
Output / Gain (2cc) Slim Tube	129 / 69



Easy Line Charger Combi BTE 2 with optional Power Pack
Dry, protect, charge!



HANSATON Charger
Small but powerful!



Rechargeable Runtime (charging from 0 % to 100 % in less than 3 hours)
Up to 24 hours: Without streaming
Up to 16 hours: AirStream technology (TV Connector/PartnerMic)
Up to 11 hours: Bluetooth® phone and music streaming

Hearing aids are shown at actual size.



HANSATON is the exclusive brand for you and your clients since 1957.

We focus on personal relationships and your success.

Our German heritage, combined with the innovative power of Sonova, enables us to create beautifully designed products, based on cutting-edge, proven technology.



hansaton.com/pro



We drive the transition to a resource-efficient and low-carbon future of our operations throughout the lifecycle of our products and services.

sonova.com/sustainability