

Hansaton Essenz

Hear the essence of life

Right from the start, your clients benefit from the advanced capabilities of HearIntelligence 6.0, with Hansaton Essenz sound processing powered by the innovative ERA™ chip.

Highlights:

- The **AI-trained** automatic program, **Intelligent AutoSurround**, ensures rapid, seamless optimization for each listening situation.
- **Bringing conversations to life in noisy environments** with StereoBeam, now including Intelligent NoiseReduction and Directional SpeechBoost.
- **Let quiet voices be heard** with Soft SpeechBoost, reducing listening effort, minimizing fatigue, and improving speech intelligibility for distant and soft talkers.¹



sound E-R shown at actual size

Wireless rechargeable RIC, with IP68 water resistance

¹ Hansaton Research Insights: Better Hearing in Quiet (2025).

² All comparisons are to the previous generation of Hansaton products.

³ Up to 16 hours, including 8 hours everyday listening, 4 hours Bluetooth® streaming, and 4 hours TV streaming. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.

Universal connectivity:

- **Stable connectivity** to enjoy uninterrupted hands-free calls and streaming.
- **Ultra-responsive ERA chip** to allow effortless switching between your client's favorite music and their favorite person.
- **Extended range** for connectivity up to 2x the distance enabling a whole new freedom of movement.²







Energy for the whole day!³



Easy Line Charger / ChargerGo RIC S
with magnetic hearing aid insertion

 **HANSATON**

Performance Levels

	9	7	5	3	1
UPDATED HearIntelligence 6.0 with Hansaton Essenz sound processing	✓	✓	✓	✓	✓
 Intelligent AutoSurround					
Intelligent AutoSurround	8	6	4	2	
Conversation in Loud Noise	✓				
Conversation in a Crowd	✓				
Conversation in a Small Group	✓	✓			
Music	✓	✓			
Quiet	✓	✓	✓		
Noise	✓	✓	✓		
Conversation in Noise	✓	✓	✓	✓	
Conversation in Quiet	✓	✓	✓	✓	✓ ¹
 Focusing on Speech					
UPDATED StereoBeam	✓	✓ ²			
Dynamic SpeechBeam	✓				
SpeechBeam		✓			
Dynamic SpeechBeam in car	✓				
UPDATED Directional SpeechBoost	✓				
UPDATED Soft SpeechBoost	✓	✓	✓		
UPDATED SpeechBoost	✓	✓	✓	✓	✓
 Localization of Acoustic Signals					
SphereSound Natural	✓				
SphereSound	✓	✓	✓		
 Comfort in Noise					
Intelligent NoiseReduction	✓	✓			
UPDATED NoiseReduction	✓	✓	✓	✓	✓
UPDATED Active Wind Block			✓		
 Personalized Adjustments					
Channels (G/AGC)	20	20	16	12	12
UPDATED Lifestyle Analyzer			✓		
NEW AI Dome Proposer			✓		
 Wireless Connectivity					
BiPhone	✓	✓	✓	✓	
Tap control	✓	✓	✓	✓	
UPDATED Universal Connectivity			✓		
UPDATED Handsfree Bluetooth® phone calls			✓		
UPDATED Bluetooth® media streaming			✓		
UPDATED stream remote App			✓		
NEW Find my Hearing Aids			✓		

Available in all performance levels

Intelligent AutoSurround: AutoStream (SpeechStream, MusicStream)

Focusing on Speech: SoundRestore frequency compression

Localization of Acoustic Signals: Pinna Effect

Comfort in Noise: Sound Impulse Manager, FeedbackManager, Tinnitus Manager

Personalized Adjustments: FocussedFit, Intelligent Acclimatization, DataLogging

Wireless Connectivity: BiLink, RogerDirect™. Optional TV Connector, PartnerMic, and Remote Control accessories

