



Focus on what's important to you!

Whether you're a music lover, an outdoor enthusiast, a social person, or a busy professional, our range of performance levels has got you covered. With HANSATON FOKUS, we offer **multiple performance levels**, so you can find the right match for your **unique listening lifestyle**.

HANSATON FOKUS hearing aids provide **great sound performance**. Plus, hearing aid options with **rechargeable Lithium-ion batteries** and **wireless connectivity**, in designs you will love.

Selecting the right performance level is crucial to maximizing your hearing experience. With our wide range of performance levels, you can rest assured that you will find **the right level that fits your unique needs**.

Don't let hearing loss hold you back from
enjoying the moments that matter.

Be in the moment...
with HANSATON FOKUS



AQ sound FS R shown at actual size.

HANSATON FOKUS Performance Levels



HANSATON is the exclusive brand for you and your clients since 1957. We focus on personal relationships and your success. Our German heritage, combined with the innovative power of Sonova, enables us to create beautifully designed products, based on cutting-edge, proven technology.



hansaton.com



We drive the transition to a resource-efficient and low-carbon future of our operations throughout the lifecycle of our products and services.

sonova.com/sustainability



Sonova AG · Laubisrütistrasse 28
CH 8712 Stäfa · Switzerland



028-6729-02/Y1.00/2023-02/MH/©2023 Sonova AG. All rights reserved.

HANSATON FOKUS PERFORMANCE LEVELS



A Sonova brand

The right level for every need

Each performance level includes the features from the lower level plus additional features.

HEAR Intelligence™

Be in the moment... with HANSATON FOKUS

Find the right level to help you communicate with confidence with more advanced HearIntelligence™ features at each performance level.

9	<div>Parties</div> <div>Conversations in the car</div> <div>Conversations in crowds or loud places</div> <div>Sporting events</div> <div>Travel</div>	LEVEL 9: FIRST Spotlight the speaker from any direction. Speech is highlighted further than any other level. Specific feature for conversations in the car from all directions around you.	Conversation in a crowd Conversation in loud noise Dynamic SpeechBeam with Intelligent NoiseReduction StereoBeam Dynamic SpeechBeam in car SphereSound Natural ConversationOptimizer Additional NoiseReduction and ConversationLift
7	<div>Restaurants</div> <div>Business meetings</div> <div>Music – live and recorded</div> <div>Shopping</div>	LEVEL 7: BUSINESS Automatically detects where speech is coming from for focus to the front or sides. Additional features reduce noise beside and behind. Music is treated specially.	Conversation in a small group Music environment StereoBeam in a manual program SpeechBeam with Intelligent NoiseReduction Additional NoiseReduction and ConversationLift
5	<div>Cafes</div> <div>Noisy transportation</div> <div>Most TV programs</div> <div>Conversation partners with quiet voices</div>	LEVEL 5: COMFORT Automatically helps with conversations in some noise. The left and right hearing aids work together to recreate a more natural listening experience. Boosts quiet voices.	Noise environment Quiet environment SphereSound ConversationBoost Additional NoiseReduction and ConversationLift
3	<div>Small group dinners in quiet</div> <div>Landline phone calls</div> <div>Games with a few people</div>	LEVEL 3: ECONOMY The hearing aids automatically adjust to loud and calm situations plus an additional feature that assists with landline phone calls.	Conversation in noise Conversation in quiet BiPhone Tap control for rechargeable hearing aids
1	<div>1 to 1 conversations</div> <div>Wireless connectivity</div> <div>Rechargeable options</div>	LEVEL 1: LOWER ECONOMY Basic functionality plus wireless connectivity and rechargeable options.	NoiseReduction ConversationLift Sound Impulse Manager Active Wind Block

