

# BE IN THE MOMENT

WITH HANSATON FOKUS  
RIC HEARING AIDS



A Sonova brand

 **HANSATON**  
*hearing & emotions*

# HANSATON FOKUS RIC RECHARGEABLE HEARING SYSTEMS

## Small sizes, great technology, stylish designs



HANSATON FOKUS RIC hearing aids come with a modern, unique design and the exceptional FOKUS sound performance. They're available in five performance levels, which include up to eight automatic listening environments, part of the innovative HearIntelligence™ technology. They also offer connectivity to all mobile phones\*, hands-free phone calls, and the freedom of rechargeability.



Sterling  
Silver



Metallic  
Anthracite




Carbon  
Black



Cashmere  
Beige



Mocha  
Brown



### Receiver type

Output / Gain (2cc)

**S**

111 / 47

**M**

114 / 51

**P**

122 / 59

**UP**

132 / 71



"The great sound performance and connectivity you want, in a design you will love. Plus easy handling, fast charging, and less waste."

CLIENT BENEFIT

## HANSATON CHARGING SYSTEMS

# More freedom: Up to 24 hours of runtime

At the heart of HANSATON FOKUS rechargeable hearing systems is a proven Lithium-ion battery with up to 24 hours of runtime and only three hours of charging time.



### Comfort Charger Combi

This Comfort Charger Combi is a charger and protective hard case in one. With the help of optional HADEO drying capsules, it also dries hearing systems when closed. With the optional Comfort Power Pack it can charge hearing systems up to seven times.

**Dry, protect, charge!**



### HANSATON Charger

The smaller HANSATON Charger combines charging comfort and attractive design. With a maximum charging time of three hours, hearing systems with Lithium-ion batteries can be charged.

**Small but powerful!**



### HANSATON T Charger

With the new and compact HANSATON T Charger, hearing systems with Lithium-ion batteries can be charged within three hours. Thanks to magnetic slots, the hearing aids can be inserted and removed effortlessly.

**Small and easy!**




### Rechargeable Runtime (charging from 0 % to 100 % in less than 3 hours)

Up to 24 hours: Without streaming

Up to 16 hours: AirStream technology (TV Connector / PartnerMic)

Up to 11 hours: Bluetooth® phone and music streaming

\* Smartphones and traditional mobile phones with compatible Bluetooth® profiles.  
The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.



HANSATON is the exclusive brand for you and your clients since 1957.

We focus on personal relationships and your success.

Our German heritage, combined with the innovative power of Sonova, enables us to create beautifully designed products, based on cutting-edge, proven technology.



[hansaton.com/pro](https://hansaton.com/pro)



We drive the transition to a resource-efficient and low-carbon future of our operations throughout the lifecycle of our products and services.

[sonova.com/sustainability](https://sonova.com/sustainability)