

Style RIC rechargeable hearing systems with HANSATON FOKUS technology



HANSATON is the exclusive brand for you and your clients since 1957.

We focus on personal relationships and your success.

Our German heritage, combined with the innovative power of Sonova, enables us to create beautifully designed products, based on cutting-edge, proven technology.




hansaton.com/pro



We drive the transition to a resource-efficient and low-carbon future of our operations throughout the lifecycle of our products and services.

sonova.com/sustainability



 Sonova AG · Laubisrütistrasse 28
CH 8712 Stäfa · Switzerland



028-6759-02/V1.00/2023-06/MH/©2023 Sonova AG. All rights reserved.

BE IN THE MOMENT

WITH HANSATON FOKUS STYLE RIC HEARING AIDS



GERMAN
DESIGN
AWARD
SPECIAL
2023

A Sonova brand





AWARD-WINNING DESIGN

Freedom meets style.

The beautiful, slender Style RICs have received the German Design Award. Their unique left and right receiver angles were designed for comfort and retention.

HIGH PERFORMANCE

based on HANSATON FOKUS technology

EASY CONNECTIVITY

to all mobile phones and tablets*

RELIABLE RECHARGEABILITY

lasting all day on a single charge in a stylish charger

HANSATON S Charger Designed to impress

Style RICs are delivered with the **elegant and easy-to-use HANSATON S Charger**, which is unique to the Style RIC. **Magnetic hold** for the hearing aids provides **easier insertion and removal** and the Style RICs are **charged in less than 3 hours**.



HANSATON FOKUS STYLE RIC RECHARGEABLE HEARING SYSTEMS

Performance meets style.




Based on the HANSATON FOKUS platform: with the outstanding performance and modern shape of the Style RIC, your clients can focus on what's most important, supported by direct connectivity, including tap control. With key features such as StereoBeam and Dynamic SpeechBeam in car they can lead an active lifestyle and enjoy social activities.

The award-winning Style RICs are available at performance levels 9, 7, 5 and 3.



AQ sound FS S

The beautiful, exclusive, slender design of the AQ sound FS S has a push-button and on-board LEDs for easy understanding of the charging status of each hearing aid. It measures 3.9 cm and is powered by a Lithium-ion battery, that lasts all day.



Receiver type	S	M	P
Output / Gain (2cc)	111 / 47	114 / 51	122 / 59

* Smartphones and traditional mobile phones with compatible Bluetooth® profiles. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. Products are shown at actual size.