

# UV Dry Box

## Electronic Cleaning and Drying Station for Hearing Aids

### Why is cleaning and drying your hearing aids important?

Maintaining your hearing systems is important, as hearing systems comprise many tiny components and highly sophisticated electronics. To safeguard the functionality of your hearing systems and to prevent damage, it is important to know how to clean and care for your hearing systems correctly. By doing so you can significantly increase the functionality and lifespan of your hearing aids. Care that quickly pays for itself.

The Clean and Dry Box cleans the hearing aid by means of UVC and dries by means of heating and airflow. The UVC light eliminates dirt and the presence of up to 99.9% of bacteria from your hearing aids. This allows the hearing aids to perform better, so that it creates a better hearing effect and can help to prevent skin irritations, itching, rashes or ear infections.



Suitable for rechargeable hearing aids.



# Cleaning & Drying

UV DRY BOX

up to 99.9%  
effective



For additional information regarding HADEO  
as well as HANSATON hearing systems and services,  
please visit us at [www.hansaton.com](http://www.hansaton.com)



 **Legal Manufacturer:**  
Wake Up International B.v.  
Bergweg 49, 3707 Aa Zeist  
The Netherlands | +31(0)30697 32 40

**Distribution:**  
Sonova AG  
Laubisruetistrasse 28  
8712 Staefa · Switzerland



# Advantages of the UV Dry Box for hearing aids

- The UVC light eliminates dirt and the possible presence of bacteria from your hearing aids, up to 99.9% effective
- Can help to prevent skin irritations, itching, rashes or ear infections
- Suitable for rechargeable hearing systems
- Box designed for easy grip and handling
- UVC is environmentally friendly, clean and efficient
- Reliable cleaning and drying
- Easy to use
- Can extend the life of hearing aids

## Cleaning and Drying Hearing Aids with HADEO UV Dry Box

The Dry Box gently dries hearing aids without chemical additives and uses very little energy. Additionally, the UVC light eliminates germs by up to 99.9%.

Removable basket holds the hearing aids with one area for the left hearing aid and one for the right

- ✓ Easy handling of the Dry Box and your hearing aids

Air flows at 30 degrees Celsius along the sides of the box for reliable drying

- ✓ Compatible with rechargeable hearing aids

5 LED lights indicate the remaining time in the cycle and the box turns off automatically after completion

- ✓ Clear and easy use



UVC light for cleansing and the elimination of germs up to 99.9%

- ✓ Can help reduce skin irritations

One-button power on/off



One program



UVC



Heat function



USB-C



Status display