

# HearIntelligence™ supports natural hearing

HearIntelligence™ technology was developed by HANSATON to support natural hearing by enhancing natural functions and processes in human hearing.

SHD 3	
<b>AutoSurround</b>	3 environments
Conversation in noise	✓
Quiet	✓
Conversation in quiet	✓
SurroundSupervisor	✓
<b>Focusing on Speech</b>	
Adaptive Directional*	✓
Fixed Directional*	✓
Pinna Effect*	✓
SurroundOptimizer	✓
ConversationLift	✓
SoundRestore	✓
<b>Comfort in Noise</b>	
NoiseReduction	✓
Sound Impulse Manager	✓
Active Wind Block	✓
FeedbackManager	✓
Tinnitus Manager	✓
<b>Personalized Adjustments</b>	
Channels (Gain / AGC)	8
FocussedFit	✓
AcclimatizationManager	✓
DataLogging	✓
<b>Wireless Connectivity</b>	
BiLink*	✓
BiPhone*	✓
PhoneConnect	✓



The housings of SHD 3 BTE and RIC hearing systems – with their special plasma coating – are dust and water resistant.

\* Not available with the jam SHD 3-10 ITE

Technical and audiological information  
SHD 3 hearing systems



HANSATON is the exclusive brand for you and your clients since 1957. We focus on personal relationships and your success. Our German heritage, combined with the innovative power of Sonova, enables us to create beautifully designed products, based on cutting-edge, proven technology.



[hansaton.com/pro](https://hansaton.com/pro)



We drive the transition to a resource-efficient and low-carbon future of our operations throughout the lifecycle of our products and services.

[sonova.com/sustainability](https://sonova.com/sustainability)



Sonova AG · Laubisrütistrasse 28  
CH 8712 Stäfa · Switzerland



028-6785-02 / V1.00 / 2023-08-08 / ©2023 Sonova AG. All rights reserved.

# SHD 3

ACOUSTICALLY IMPRESSIVE



# SHD 3 hearing aids are aesthetically appealing and acoustically impressive

## Boosts conversations and suppresses noise.

People with normal hearing can ignore disturbing noise and concentrate instead on the important sounds they want to hear. HearIntelligence™ technology supports your clients' hearing by processing incoming sounds and reducing them to essential components that are either boosted or suppressed, providing an optimized hearing experience.

## Design meets performance.

HANSATON combines high performance with stylish design in each of the SHD 3 styles. You can select from a RIC design that is very high-fashion, or BTEs that are known for their reliability, or ITEs that can be custom tailored to your clients' needs. Stylish designs can make your clients feel comfortable trying hearing aids. And the housings of the SHD 3 BTE and RIC hearing systems – with their special plasma coating – are dust and water resistant, making them even more flexible to match your clients' lifestyles.

## Sound Optimization.

- Based on proven technology, SHD 3 covers a wide range of hearing losses and has appealing designs to meet your clients' needs
- With the FocussedFit workflow you can get your clients on their way to hearing better quickly
- The AcclimatizationManager starts them at a level that is comfortable and provides a way to reach an optimal setting over time as they become accustomed to amplification

Help your clients hear the important moments in their lives.

Small, custom-made and discreet.

## ITE hearing aids

SHD 3 ITE hearing systems deliver great quality and modern technology combined in small custom designs.


Optional functions:  
Push button, Volume control,  
Multi-microphone technology, Telecoil



Beautifully designed and compact.

## RIC hearing aid


The SHD 3 RIC hearing system combines top quality and innovative technology with appealing, discreet design. A compelling option.



Reliable and robust.


## BTE hearing aids

The SHD 3 BTE hearing systems provide a winning combination of great technology, robust housings and convenient handling.




Receiver Types	Output / Gain (2cc)	jam SHD 3-10	jam SHD 3-312 Dir W
Moderate	109/40	(✓)	(✓)
Power	115/50	✓	✓
Super Power	119/60	(✓)	(✓)
Ultra Power	127/70		(✓)
Shell styles		IIC/CIC/MC/ITC/HS/FS	ITC/HS/FS

✓ Default (standard)  
(✓) Optional



Receiver 3.0	Standard (xS)	Power (xP)	Super Power (xSP)
Output / Gain (2cc)	111 / 47	124/57	125/62
Open Dome	✓	✓	
Closed Dome	✓	✓	
Power Dome	✓	✓	
SlimTip	✓	✓	
cShell	✓	✓	✓



Receiver 3.0	jam SHD 3-RS13	beat SHD 3-RS675
Output / Gain (2cc) Earhook	131/66	
Output / Gain (2cc) Slim Tube	126/58	
Output / Gain (2cc) Filtered Earhook		133/75
Output / Gain (2cc) Unfiltered Earhook		139/82
Output / Gain (2cc) Power Slim Tube		139/82

All products shown at actual size