

HANSATON ITE product portfolio October 2017.

NEW
with SHD
Technology

NEW



NEW



NEW




NEW



NEW



NEW



IIC (10 Omni)

CIC (10 Omni/Omni W)

Mini Canal (10 Omni W/Dir W)

Canal (312 Dir W)

Half Concha (312 Dir W/13 Dir W)

Concha (312 Dir W/13 Dir W)

Names	Performance Level				Battery	Power	Accessories	Colors
jam SHD IIC (10 Omni)	-	7	5	3	10	M, P, SP	IIC remote	-
jam SHD CIC (10 Omni/Omni W)	-	7	5	3	10	M, P, SP	RCV2*, uDirect3*, uStream*, uTV3*, uMic2*	-
jam SHD Mini Canal (10 Omni W/Dir W)	9	7	5	3	10	M, P, SP	RCV2, uDirect3, uStream, uTV3, uMic2	-
jam SHD Canal (312 Dir W)	9	7	5	3	312	M, P, SP, UP	RCV2, uDirect3, uStream, uTV3, uMic2	-
jam SHD Half Concha (312 Dir W/13 Dir W)	9	7	5	3	312 / 13	M, P, SP, UP	RCV2, uDirect3, uStream, uTV3, uMic2	-
jam SHD Concha (312 Dir W/13 Dir W)	9	7	5	3	13 / 312	M, P, SP, UP	RCV2, uDirect3, uStream, uTV3, uMic2	-

Names	Performance Level	L0max/Gmax (2ccm-coupler)	Channels (G/AGC)	No. of programs (AutoSurround SHD/Manual/Wireless)	Battery Size	AutoSurround SHD	AutoSurround SHD Quiet	AutoSurround SHD Conversations in quiet	AutoSurround SHD Conversations in noise	AutoSurround SHD Conversations in a crowd	AutoSurround SHD Conversations in a small group	AutoSurround SHD Noise	AutoSurround SHD Music	SurroundSupervisor SHD	SphereSound SHD Dynamic	SphereSound SHD Personalized	SphereSound SHD Static	SpeechBeam SHD-3	SpeechBeam SHD-2	SpeechBeam SHD-1	SurroundOptimizer SHD	AcclimatizationManager	ConversationLift+ / SpeechLift (SL)	NoiseReduction	FeedbackManager	Sound Impulse Manager SHD	Active Wind Block	SoundRestore	DataLogging	Tinnitus Manager	BiLink	BiPhone	Bluetooth (accessory needed)	PhoneConnect	Pushbutton optional	Telecoil optional	
IIC (10 Omni)	7	M: 109/40, P: 115/50, SP: 119/60	16 7/3/3	10																																	
IIC (10 Omni)	5	M: 109/40, P: 115/50, SP: 119/60	12 6/3/3	10																																	
IIC (10 Omni)	3	M: 109/40, P: 115/50, SP: 119/60	8 3/3/3	10																																	
CIC (10 Omni/Omni W)	7	M: 109/40, P: 115/50, SP: 119/60	16 7/3/3	10																													*	*			
CIC (10 Omni/Omni W)	5	M: 109/40, P: 115/50, SP: 119/60	12 6/3/3	10																													*	*			
CIC (10 Omni/Omni W)	3	M: 109/40, P: 115/50, SP: 119/60	8 3/3/3	10																													*	*			
Mini Canal (10 Dir W)	9	M: 109/40, P: 115/50, SP: 119/60	20 7/3/3	10																																	
Mini Canal (10 Omni W/Dir W)	7	M: 109/40, P: 115/50, SP: 119/60	16 7/3/3	10																																	
Mini Canal (10 Omni W/Dir W)	5	M: 109/40, P: 115/50, SP: 119/60	12 6/3/3	10																																	
Mini Canal (10 Omni W/Dir W)	3	M: 109/40, P: 115/50, SP: 119/60	8 3/3/3	10																																	
Canal (312 Dir W)	9	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	20 7/3/3	312																																	
Canal (312 Dir W)	7	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	16 7/3/3	312																																	
Canal (312 Dir W)	5	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	12 6/3/3	312																																	
Canal (312 Dir W)	3	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	8 3/3/3	312																																	
Half Concha (312 Dir W/13 Dir W)	9	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	20 7/3/3	312																																	
Half Concha (312 Dir W/13 Dir W)	7	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	16 7/3/3	312																																	
Half Concha (312 Dir W/13 Dir W)	5	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	12 6/3/3	312																																	
Half Concha (312 Dir W/13 Dir W)	3	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	8 3/3/3	312																																	
Concha (312 Dir W/13 Dir W)	9	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	20 7/3/3	13																																	
Concha (312 Dir W/13 Dir W)	7	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	16 7/3/3	13																																	
Concha (312 Dir W/13 Dir W)	5	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	12 6/3/3	13																																	
Concha (312 Dir W/13 Dir W)	3	M: 109/40, P: 115/50, SP: 119/60, UP: 127/70	8 3/3/3	13																																	

* Only for Omni W