

HANSATON BTE product portfolio May 2018.

NEW



AQ jam SHD RS

jam SHD S312

jam SHD RS13

jam SHD S13

Names	Performance Level				Battery	Power	Accessories	Colors
AQ jam SHD RS (Rechargeable BTE-System with pushbutton)	9	7	5	3	Li-Ion	P	RCV2, uDirect3, uStream, uTV3, uMic2	1 – 8
jam SHD S312 (Slimtube BTE-System with pushbutton)	9	7	5	3	312	M	RCV2, uDirect3, uStream, uTV3, uMic2	1 – 8
jam SHD RS13 (BTE-System with pushbutton)	9	7	5	3	13	P	RCV2, uDirect3, uStream, uTV3, uMic2	1 – 7
jam SHD S13 (BTE-System with pushbutton)	9	7	5	3	13	P	RCV2, uDirect3, uStream, uTV3, uMic2	1 – 7



(1) White Pearl



(2) Cloudy Grey



(3) Sterling Silver



(4) Space Titan



(5) Carbon Black



(6) Sparkling Bronze



(7) Sandy Beige



(8) Black Edition
for Performance Level 9 and 7

HEAR IntelligenceTM

Names	Performance Level	L0max/Gmax (2ccm-coupler)	Channels (G/AGC)	No. of programs (AutoSurround SHD/Manual/Wireless)	Battery Size	AutoSurround SHD	AutoSurround SHD Quiet	AutoSurround SHD Conversations in quiet	AutoSurround SHD Conversations in noise	AutoSurround SHD Conversations in a crowd	AutoSurround SHD Conversations in a small group	AutoSurround SHD Noise	AutoSurround SHD Music	SurroundSupervisor SHD	SphereSound SHD Dynamic	SphereSound SHD Personalized	SphereSound SHD Static	SpeechBeam SHD-3	SpeechBeam SHD-2	SpeechBeam SHD-1	SurroundOptimizer SHD	AcclimatizationManager	ConversationLift+ / SpeechLift (SL)	NoiseReduction	FeedbackManager	Sound Impulse Manager SHD	Active Wind Block	SoundRestore	DataLogging	Tinnitus Manager	BiLink	BiPhone	Bluetooth (accessory needed)	PhoneConnect	Audio input	Pushbutton	Telecoil	IP68		
AQ jam SHD RS	9	Earhook: 131/66 Slimtube: 126/58	20	7/3/3	Li-Ion	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S312	9	Earhook: 129/63 Slimtube: 122/56	20	7/3/3	312	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD RS13	9	Earhook: 131/66 Slimtube: 126/58	20	7/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S13	9	Earhook: 131/63 Slimtube: 123/55	20	7/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
AQ jam SHD RS	7	Earhook: 131/66 Slimtube: 126/58	16	7/3/3	Li-Ion	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S312	7	Earhook: 129/63 Slimtube: 122/56	16	7/3/3	312	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD RS13	7	Earhook: 131/66 Slimtube: 126/58	16	7/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S13	7	Earhook: 131/63 Slimtube: 123/55	16	7/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
AQ jam SHD RS	5	Earhook: 131/66 Slimtube: 126/58	12	6/3/3	Li-Ion	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S312	5	Earhook: 129/63 Slimtube: 122/56	12	6/3/3	312	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD RS13	5	Earhook: 131/66 Slimtube: 126/58	12	6/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S13	5	Earhook: 131/63 Slimtube: 123/55	12	6/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
AQ jam SHD RS	3	Earhook: 131/66 Slimtube: 126/58	8	3/3/3	Li-Ion	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S312	3	Earhook: 129/63 Slimtube: 122/56	8	3/3/3	312	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD RS13	3	Earhook: 131/66 Slimtube: 126/58	8	3/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.
jam SHD S13	3	Earhook: 131/63 Slimtube: 123/55	8	3/3/3	13	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.	.